



Compact, portable and exceptionally versatile



The Oxford® Journey lift is an ultra compact, portable folding stand aid with best-in-class functionality. It is compact and lightweight allowing care-givers to manoeuvre in tight spaces and can be folded for easy storage or transportation.

The unique adjustable cow-horn mechanism supports a wider range of patient heights and sizes, allowing your facility to have one lift that meets the needs of many.



# Key Features

- Safe Working Load: 155kg/341lbs/24st
- Folding Design
- Smart™ Monitor Controls
- Ergonomic Design Focus
- Unique Adjustable Cow-horn
- Height Adjustable Knee Pad
- Removable Foot Tray





# Folding Design

Unlike stand aids on the market, the Journey can fold making it easy to transport and store.

This functionality gives active users the option to travel outside the confines of their own home.



# Smart™ Monitor Controls

The Smart™ Monitor from Oxford is a state-of-the-art control system providing one of the industry's most accurate lift counters.

The intelligent system advises the user to schedule key maintenance routines, ensuring your lift remains in perfect working order.

Additionally, the Smart™ Monitor alerts the user and service engineer to potential misuse (attempted overloads above safe working load).



# Ergonomic Design Focus

The Journey was designed with the 5 key principles of moving and handling in mind. This ensures Oxford's lift designs are simple, safe and above all comfortable to use.

1. Get close to the load
2. Use a wide stable base
3. Ensure a comfortable firm grasp
4. Keep spine close to neutral
5. Make sure movement is smooth





# Unique Adjustable Cow-horn

The Journey's unique adjustable cow-horn ensures maximum flexibility and improved patient comfort.

With three cow-horn height options, the Journey can support a wider range of patient heights and sizes.





# Height Adjustable Knee Pad

The Journey's contoured knee pad is height adjustable to ensure the patient feels comfortable and fully supported during transfer.

A retractable safety belt provides additional support for patients needing further reassurance.



# Removable Foot Tray

By removing the foot tray the Journey transforms into a rehabilitation aid that allows a patient to fully stand.

Once standing, the patient can be encouraged to walk independently.



# Sling Compatibility

The Journey can be used in conjunction with either a **Deluxe Standing** or **Deluxe Transport** sling. This provides the carer with the option of performing various care duties, including:

- ✓ Toileting
- ✓ Standing
- ✓ Point to point transfers
- ✓ Walking and rehabilitation activity

